

REPLACING THE SAFE-T-RAIL FOAM GRIP



1. Carefully cut off the old foam grip with a box cutter or knife making sure you don't chip the paint on the Safe-T-Rail.
2. Pop off the end cap of the Safe-T-Rail with a flat head screw driver.
3. Liberally spray the Safe-T-Rail with a non-soapy solution like Simple Green or 409. You can also dunk the foam grip in this solution prior to putting it on the Safe-T-Rail.
4. Pull the new foam grip on to the Safe-T-Rail until you have it in the correct position.
5. Put the end cap back on the Safe-T-Rail.